

SQUID Characterization RUN 18 and RUN 19, measurements taken on 15.05.2024 and 24.05.2024, respectively

Nicolas Sierra Manrique

Abstract—The following paper presents the results of the latest characterization of SQUID sensors for the Nucleus experiment.

Contents

1	Introduction	1
2	RUN 15	1
2.1	QUAD 190	2
3	RUN 19	2
4	Noise Code Run 18	3
5	Noise Code Run 19	4
	References	4

1. Introduction

The previous characterization indicates that the SQUID-noise patterns depend on the circuit setup. Currently, there are four configurations installed in the SQUID system:

- A SQUID with a short installed (see Figure 1).
- A SQUID connected to the breakout box but not to a detector (see Figure 2), also referred to a Shunt-shortened circuit.
- A SQUID with a bias current flowing entirely through the shunt resistances (see Figure 3), also referred as an Open circuit.
- A SQUID with a TES detector in operation (see Figure 4).

		SQUID				
		1	2	3	4	
		5	6	7	8	
		Breakout box				
SCSI 1	QUAD 185	O	O	O	O	
	QUAD 186	O	SH	SH	SH	TUM
SCSI 3	QUAD 190	SH	SH	D	SH	STAR
	QUAD 191	SH	SH	SH	SH	STAR
SCSI 2	QUAD 107	S	S	S	S	STAR
	QUAD 189	S	S	S	S	STAR

Table 1. Measured detector configuration **O** stands for Open, **SH** stands for Shunt, **S** stands for Short, and **D** for Detector




Figure 1. Shorted Circuit




Figure 2. Circuit connected to the Shunt resistances, without a detector








Figure 3. Open Circuit, entire Bias current flowing through the shunt resistances

2. RUN 15

The SQUID noise of QUADS 107 and 189 with STAR Cryogenic wiring was measured with the Pulse Tube on and with the Pulse Tube off.

Table 2. Amplitude of the SQUID patterns		
	SQUID 1 [V]	SQUID 2 [V]
QUAD 107	4,680	5,893
QUAD 189	3,459	4,184
	SQUID 3 [V]	SQUID 4 [V]
QUAD 107	3,591	3,545
QUAD 189	4,870	4,372

Both QUADS (107, 189) were installed with a short attached to the breakout box

**Figure 4.** Circuit with a TES Detector**Figure 5.** Noise of QUAD 107 with pulse tube ON**Figure 6.** Noise of QUAD 189 with pulse tube ON**Figure 7.** Noise of QUAD 189 and QUAD 107 with pulse tube OFF**Table 3.** Mean SQUID noise for QUAD 107 and 189 around 10 Hz, with pulse tube ON


	SQUID 1 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 2 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 107	118	84
QUAD 189	150	125
	SQUID 3 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 4 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 107	6	114
QUAD 189	146	97

Table 4. Mean SQUID noise for QUAD 107 and 189 around 10 Hz, with pulse tube OFF

	SQUID 1 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 2 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 107	126	95
QUAD 189	N/A	N/A
	SQUID 3 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 4 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 107	5	N/A
QUAD 189	N/A	117

2.1. QUAD 190

Additionally, the noise of SQUID 3 from QUAD 190 with pulste on and off was measured.

**Figure 8.** Noise of SQUID 3 from QUAD 190 with pulste on and off**Table 5.** Mean SQUID noise of SQUID 3 in QUAD 190 around 10 Hz, with pulse tube ON and OFF

	SQUID 3 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
ON	3.8
OFF	3.7

3. RUN 19

In Run 19 the noise of QUADs 185, 186, 190 and 191 was measured.




Figure 9. Noise of QUAD 185




Figure 12. Noise of QUAD 191

Table 6. Mean SQUID noise at 10 Hz for QUAD 185, 186, 190 and 191 from RUN 19

	SQUID 1 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 2 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 185	275	82
QUAD 186	150	9
QUAD 190	613	451
QUAD 191	26	75

	SQUID 3 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 4 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 185	10	31
QUAD 186	57	64
QUAD 190	3.7	2.7
QUAD 191	95	45




Figure 10. Noise of QUAD 186

Table 7. Mean SQUID noise for each QUAD around 10 Hz, from RUN 8 cells in blue are the SQUIDs in which a detector is connected: Table from [1]

	SQUID 1 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 2 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 185	133.50	8.6462
QUAD 186	125.73	15.381
QUAD 190	150.22	262.83
QUAD 191	4.8968	4.1133

	SQUID 3 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$	SQUID 4 $\frac{\text{pA}}{\sqrt{\text{Hz}}}$
QUAD 185	7.0264	98.251
QUAD 186	100.79	6.7787
QUAD 190	5.1965	7.4317
QUAD 191	8.9389	5.9237




Figure 11. Noise of QUAD 190

4. Noise Code Run 18

Since the data taken from Run 18 was saved as int32, some modifications had to be made to the conversion factors to account for the VDAQ values.

```
1 nhead = 6
2 toVolt = 39.3216/((2**24) - 1)
3 touAmp= (77/55)*toVolt
```

Code 1. ToVolt conversion value for Int32.

Here is an example code for the Fourier transform which allows for the measurement of the noise:

```
1 for i in range(4):
2     file_path = r"C:\Users\nicos\Downloads\
    clusterCryostatRun18\noise\noise_SQ"+ str(i + 1) +
    "_000.bin"
```

```

3     #print("squidchar_noise_186" + str(i+1)+"M")
4     fs=100e3
5     nperseg=int(0.5*fs)
6
7     if os.path.isfile(file_path):
8         #readRDT(file_path)
9         fd = open(file_path, 'rb')
10        read_data = np.fromfile(file=fd, dtype=np.int32
11        )
12        fd.close()
13        read_data= toVolt*read_data[12:]/1.70
14        f, pxx = sig.welch(read_data, fs=fs, nperseg=
15        nperseg)
16        f = f[:-1]
17        ppx = 1e6 * np.sqrt(ppx[:-1])
18        P10 = np.mean(ppx[(f >= 7) & (f <= 11)])
19        P103 = np.mean(ppx[(f >= 900) & (f <= 1100)])
20        print(P10)
21        plt.plot(f, ppx,label=f"SQUID_"+ str(i + 1))
22
23    else:
24        print("No data taken")
25    #plt.title("QUAD 185",fontsize= 16)
26    #plt.axhline(y=4, color='r', linestyle='--', linewidth
27    =5,label=" Shunt Noise")
28    plt.legend()
29    plt.title("QUAD 107")
30    plt.yscale("log")
31    plt.xscale("log")
32    plt.grid(True)
33    plt.xlabel(r'$f$ [Hz]', fontsize= 16)
34    plt.ylabel(r'$\frac{P}{A}\{\sqrt{Hz}\}$', fontsize= 16)
35    plt.show()

```

Code 2. Noise Code for QUAD 185.

5. Noise Code Run 19

Run 19 was saved as int16, some modifications had to be made to the conversion factors to account for the VDAQ values.

```

1 nhead = 12
2 toVolt = 39.3216/((2**16) -1)
3 touAmp= (77/55)*toVolt

```

Code 3. ToVolt conversion value for Int16.

References

- [1] N. S. Manrique, *Characterization of squid sensors for the nucleus experiment*, 2024.